



# Triangle Furniture Systems

PART OF PLANMECA GROUP



CLEAN

SC193"

SOILED

## Instrument Recirculation Procedure

BASED ON CDC, OSHA & ADA RECOMMENDATIONS

By respecting the 10-step procedure, your instruments will be clean, dry and sterile. Never forget, always wash hands before going from the soiled area (red) to the clean area (blue). **Make it a team rule!**



01  
Protective  
Equipment



02  
Waste  
Management



03  
Pre-Soaking



04  
Cleaning



05  
Rinsing



06  
Drying



07  
Lubrication



08  
Packaging



09  
Sterilization



10  
Storage



# Systemes de Mobilier Triangle

UNE DIVISION DU GROUPE PLANMECA



PROPRE

SC193"

SOUILLÉ

## Procédure de Recirculation des Instruments

BASÉE SUR LES RECOMMANDATIONS DU CDC, OSHA & ADA

En respectant la procédure en 10 points, vos instruments seront propres, secs et stériles. De plus, n'oubliez jamais qu'il faut se laver les mains avant de passer de l'environnement souillé (rouge) à l'environnement stérile (bleu). **Faites-en une règle d'équipe!**

